

Sensing Platform

Release Notes PSP_6.1.9.519-1504cf96

Rev. 1.0 — 26-Jul-2024

Accepted

Company Confidential

Document information

Info	Content
Keywords	PSP Library, Release Notes, Biosensing by PPG
Abstract	This document provides the release notes for the specified PSP Library software release

PHILIPS

Table of Contents

1 Introduction 3

1.1 Purpose..... 3

1.2 Definitions, Acronyms and Abbreviations 3

1.3 References..... 3

2 Release overview 4

2.1 Dependency overview 4

2.2 Release Content 4

3 Release PSP_ 6.1.9.519-1504cf96 for CoreKinect..... 5

3.1 Release variant 5

3.2 Change Requests 5

3.3 Problem Reports 6

3.4 API & documentation changes 6

4 Latest Release Known issues / limitations 7

4.1 Known issues (to be solved) 7

4.2 Accepted issues (NOT to be solved)..... 7

1 Introduction

1.1 Purpose

This document describes software release **PSP_6.1.9.519-1504cf96**.
It gives an overview of the implemented change requests, problem reports and known issues.

1.2 Definitions, Acronyms and Abbreviations

Abbreviation	Description
CR	Change Request, Feature
PR	Problem Report, Bug

1.3 References

Abbrev.	Title	Author(s)	Version
-	-	-	-

2 Release overview

Date	PSP version no.	Status
26 Jul 2024	PSP_6.1.9.519-1504cf96-release	6.1.9 library release. Flash Code starting address: 0x000E0000.

2.1 Dependency overview

Dependency	Identification
-	-

2.2 Release Content

Item	Description
docs	Integration documents and Release Notes.
psp.bin psp.hex	PSP library. Only one is needed.
psp.h	PSP library API header file.
fx_datatypes.h psp.c psp_map.h	FXI library API support files.

3 Release PSP_6.1.9.519-1504cf96 for CoreKinect

3.1 Release variant

This section describes the PSP library variant

Configuration	Setting
Architecture	Arm M4F
Metrics	<u>Biosensing (by PPG):</u> PSP_METRIC_ID_ACCELERATION PSP_METRIC_ID_ACTIVITYCOUNT PSP_METRIC_ID_ACTIVITYTYPE PSP_METRIC_ID_COMPRESSED_ACCELERATION PSP_METRIC_ID_COMPRESSED_PPG PSP_METRIC_ID_ENERGYEXPENDITURE PSP_METRIC_ID_FITNESSINDEX PSP_METRIC_ID_HEARTBEATS PSP_METRIC_ID_HEARTRATE PSP_METRIC_ID_LOWPOWERENERGYEXPENDITURE PSP_METRIC_ID_LOWPOWERHEARTRATE PSP_METRIC_ID_MOTIONCADENCE PSP_METRIC_ID_PPG_AMBIENT PSP_METRIC_ID_PPG_INFRARED PSP_METRIC_ID_PPG_MOTION PSP_METRIC_ID_PPG_GREEN PSP_METRIC_ID_PPG_RED PSP_METRIC_ID_PRIVATE_DATA PSP_METRIC_ID_RESPIRATIONRATE PSP_METRIC_ID_RESTINGHEARTRATE PSP_METRIC_ID_SKINPROXIMITY PSP_METRIC_ID_SLEEPSTAGES PSP_METRIC_ID_SPEED PSP_METRIC_ID_SPO2 PSP_METRIC_ID_STRESSLEVELHEARTRATE PSP_METRIC_ID_VO2MAX
Fixed Address	START_ADDRESS: 0x000E0000
Default SpO2 Calibration Coefficient	
Evaluation restrictions	Heart rate upper bond: 120 bpm SpO2 lower bond: 90%

3.2 Change Requests

This section describes the change requests that are implemented in this release.

Description	Reference	Source

3.3 Problem Reports

This section describes the problem reports that are solved in this release.

Description	Reference	Source

3.4 API & documentation changes

This section summarizes the API changes compared to the previous release.

Target	Description	Reference	Source

4 Latest Release Known issues / limitations

4.1 Known issues (to be solved)

The following issues are known and accepted in this release

Description	Reference	Source
-	-	-

4.2 Accepted issues (NOT to be solved)

The following issues will not be solved

Description	Reference	Source
-	-	-

Legal information

Disclaimers

This document may be subject to change without notice. Philips Electronics Nederland B.V. does not give any representations or warranties as to the accuracy or completeness of information included herein and shall have no liability for the consequences of use of such information.

This document is provided “as is” with no warranties, expressed or implied, including but not limited to any implied warranty of merchantability, fitness for a particular purpose, or freedom from infringement. Philips Electronics Nederland B.V. or its Affiliates may have patents or pending patent applications, or other intellectual property rights that relate to the described subject matter. Philips and its logo are trademarks of Koninklijke Philips N.V. All other names are the trademark or registered trademarks of their respective holders. The furnishing of this document does not provide any license, expressed or implied, by estoppels or otherwise, to any such patents, trademarks or copyrights, or other intellectual property rights. Philips Electronics Nederland B.V. assumes no responsibility for errors or omissions in this document; nor does Philips Electronics Nederland B.V. make any

commitment to update the information contained herein. This document is subject to change without notice.

Patents

Notice is herewith given that the subject device uses one or more of the following patents and that each of these patents may have corresponding patents in other jurisdictions.

Philips has been awarded a number of patents on the technology used in the BioSensing by PPG offering in Europe, Asia and the United States.

Trademarks

Notice: All referenced brands, product names, service names and trademarks are property of their respective owners.

CoreKinect LLC